

Advantium RX Alarm

Designed for use with the CONCOA Medical IntelliSwitch, the new Advantium Rx offers superior integration, protection, and convenience by allowing end-users to remotely monitor the status of the medical IntelliSwitch. It can serve as a single remote monitoring source or be conveniently linked to a facility alarm system or additional Advantium Rx alarms.



Advanced Features

- Specifically for use with Medical IntelliSwitch
- Five Input channels
- Five Output Channels
- Dry contact relay outputs
- Normally Open/Closed output capability
- Highly visible lighted status notification
- Distinctive Audible Notification
- Auto-reset when alarm conditions clear
- Audible Alarm Silence feature
- Ability to “Daisy Chain” multiple alarms
- Convenient on/off switch
- Integrated power supply

Specifications

- **Audio**
93db audible alarm
- **Power**
120 VAC or 220 VAC
50/60 Hz
- **Outputs**
5 channels
Dry Contacts
Normally Open/Normally Closed
(jumper selectable)
Contact Rating .5A 24V
- **Inputs**
5 channels
Accepts Normally Open/Normally Closed
contact closure (selectable)
Resettable fuse protection
- **Connectors**
Input: 8 pin circular
Output: 6 pin circular
- **Dimensions**
4 3/4”h x 7 1/2” w x 2 5/8” d

Ordering Information

Option	Order No.	Description
Remote Advantium Rx Alarm 120 VAC	574-0001	Provides visual and audible notification of the status of the Medical IntelliSwitch and Reserve Gas System
Reomote Advantium Rx Alarm 220 VAC	574-0002	Provides visual and audible notification of the status of the Medical IntelliSwitch and Reserve Gas System
Medical IntelliSwitch to Advantium Rx Cable Assembly	529-5325-XX	Pre-assembled cable assembly to connect between the Medical IntelliSwitch and the Advantium Rx Alarm. “XX” denotes cable length. Replace “XX” in part number with 03 (3’), 10 (10’), 25 (25’), 50 (50’), 100 (100’)
Reserve Gas Alarm Cable Assembly	529-5326	Pre-assembled cable assembly to connect Reserve Gas Module to Medical IntelliSwitch
Reserve Gas Alarm Extention Cable Assembly	529-5327-XX	Pre-assembled extention cable assembly to connect between the Medical IntelliSwitch and the Reserve Gas Module. Replace “XX” in part number with 25 (25’), 50 (50’), 100 (100’)
Advantium Rx Output Cable Assembly	529-5328-XX	Pre-assembled cable assembly to connect the output contacts of the Advantium Rx Alarm to other equipment. One end of cable has wires only. Replace “XX” in part number with 03 (3’), 10 (10’), 25 (25’), 50 (50’), 100 (100’)